

Do it **right** the **first** time!
Introducing a **new** solution to an old problem

DispoClean

Compatible with all Major Brands

unique
single passage cleaning process
for endoscopes



Available in two sizes:
Yellow for 3.2 - 4.2 mm Channels
Green for 2.8 - 3.2 mm Channels

DispoClean's active group of flexible blades are able to eliminate a maximum of residues.

DispoClean's multiblade technology means only a single one way pass is required to achieve reliable, faster and high quality cleaning of endoscope lumens prior to sterilisation.

EFFECTIVE

DispoClean utilises an active group of three sets of flexible multi-blades to thoroughly clean all internal endoscope channels in one simple operation.

The one piece can be used for all channels of the one scope after washing in an enzyme based cleaning solution between passes.

EFFICIENT

A single one directional pass facilitates repeatable, faster and high quality cleaning of lumens than traditional brushing.

SAFETY

DispoClean's unique proven patented process offers, safer, faster and more efficacious pre-disinfection cleaning.


COMPREHENSIVE RANGE

Available in two sizes (colour coded) for use with all endoscopes with channel diameters between 2.8mm and 4.2mm.

FORMATS

Choice of dispense formats:

- > Box of 250 individually wrapped units
- > Compact dispense reel of 250 units



CODE	DESCRIPTION	PACKAGING
DC2.8REEL	GREEN GASTRO CLEANING BLADES 2.8-3.2 CHANNEL SIZE	REEL 250
DC2.8SINGLES	GREEN GASTRO CLEANING BLADES 2.8-3.2 CHANNEL SIZE	BOX 250 INDIVIDUALLY WRAPPED
DC2.8INDIV	GREEN GASTRO CLEANING BLADES 2.8-3.2 CHANNEL SIZE	1 ONLY
DC3.2REEL	YELLOW COLON CLEANING BLADES 3.2-4.2 CHANNEL SIZE	REEL 250
DC3.2SINGLES	YELLOW COLON CLEANING BLADES 3.2-4.2 CHANNEL SIZE	BOX 250 INDIVIDUALLY WRAPPED
DC3.2INDIV	YELLOW COLON CLEANING BLADES 3.2-4.2 CHANNEL SIZE	1 ONLY

Internationally Use | Microbiologically Tested | Fast | Efficient | Effective | Safe

What will **you** do with the time **you** save?