



STERILE SURGICAL SKIN MARKER High Quality! Low Cost



The Surgical Skin Marker features a versatile tip for precise line marking, or when held at an angle a broad line marking, to allow demarcation of selected sites, prior to, or during a surgical procedure. The marker's gentian violet ink is non toxic, non tattooing and surgical prep resistant. Markings can be removed from the skin with isopropyl alcohol. Included is a 15cm, flexible ruler calibrated in both centimetres and inches. Widely used in general surgery marking, radiotherapy and B-type ultrasound positioning identification and also used in medical teaching.

Features:

- The marker's gentian violet ink is recognised by CE and has been tested by an independent laboratory
- The ink is non toxic, non tattooing and surgical prep resistant
- Medium (regular) tip
- Markings can be removed from the skin with isopropyl alcohol
- Packaged with a 15cm flexible, high quality, laminated ruler both in centimetres and inches
- Meets the single use requirement to avoid infection



SMPARMT Regular is a Medium Tip (purple lid)



Packaging Sample



Medium (Regular) Tip

CODE	DESCRIPTION
SMPARMT	LIBERTY SKIN MARKER & RULER, STERILE MEDIUM TIP (BX 25)

E&OE First Published: 8 August 2009
Last Modified: 4 January 2011
File Name: info_SMPAR_SURGICAL-SKIN-MARKER

