



HOT & COLD PACKS

INSTANT COLD COMPRESS - UREA

Usage

Simply squeeze the pouch like the instant hot compress. The contents become immediately cold and remain cold for 14-40 minutes and achieves a low temperature of -4°C

Characteristic

Produce becomes cold on its own after activation. Convenient and portable, for use whenever and wherever cold therapy is needed.

Details

Size: H-228mm x W-152mm
Weight: 260g
Code: INTICEUREA
Description: INSTANT COLD COMPRESS (BOX 40)
Pack Size: 40 per carton



REUSEABLE HOT/COLD GEL PACK IN CLOTH

Usage

To heat, microwave or boil by following instructions on back of the compress.

To obtain cold therapy, freeze for at least 1 hour.

This gel pack can reach a high of 60°C and a low of -12°C . Maintains therapy for approximately 30 minutes.

Details

Size: H-125mm x W-250mm
Code: INTGEL3510
Description: REUSABLE HOT/COLD GEL PACK IN CLOTH
Pack Size: 30 per carton



REUSEABLE HOT/COLD GEL PACK

Usage

To heat, microwave or boil by following instructions on back of the compress.

To obtain cold therapy, freeze for at least 1 hour. Can reach a high of 60°C and a low of -12°C . Maintains therapy for approximately 30 minutes.

Characteristic

Can be reused for either cold or hot. Remains soft and pliable after freezing.

Cleans easily with soap and water.

Details

Size: H-175mm x W-125mm
Weight: 110g
Code: INTGEL3057
Description: REUSABLE HOT/COLD GEL PACK
Pack Size: 100 per carton





HOT & COLD PACKS

BODY WARMER SELF ADHESIVE

Usage

Open the outer package and shake gently, leaving the internal pack to activate with the exposure to the air. It will gradually get warm and once temperature is reached it will last for around 12-15 hours, providing warmth for the user. Temperature can reach 50°C-68°C.

Characteristic

Product provides long lasting warm to assist user in cold weather. Comes with adhesive to apply to body or clothing.

Details

Size: H-175mm x W-125mm
Code: INTBODYWARMER
Description: BODY WARMER SELF ADHESIVE
Pack Size: 240 per carton



MAGIC CAT BODY WARMER

Usage

Open the outer package and shake gently, leaving the internal pack to activate with the exposure to the air. It will gradually get warm and once temperature is reached it will last for around 12-15 hours, providing warmth for the user. Temperature can reach 50°C-68°C.

Characteristic

Product provides long lasting warm to assist user in cold weather. Comes with adhesive to apply to body or clothing.

Details

Size: H-145mm x W-110mm
Code: INTMAG4046CAT
Description: INTCO MAGIC WARMER SELF ADHESIVE - CAT SHAPED
Pack Size: 50 per carton

