

## SURGICAL SKIN MARKER

The Surgical Skin Marker features a versatile tip for precise line marking, or when held at an angle a broad line marking, to allow demarcation of selected sites, prior to, or during a surgical procedure. The marker's gentian violet ink is non toxic, non tattooing and surgical prep resistant. Markings can be removed from the skin with isopropyl alcohol. Included is a 15cm, flexible ruler calibrated in both centimetres and inches. Widely used in general surgery marking, radiotherapy and B-type ultra-sonography positioning identification and also used in medical teaching.

*Liberty*<sup>®</sup>

### Features:

- The marker's gentian violet ink is recognised by CE and has been tested by an independent laboratory
- The ink is non toxic, non tattooing and surgical prep resistant
- Regular tip (0.5mm line)
- Markings can be removed from the skin with isopropyl alcohol
- Packaged with a 15cm flexible ruler both in centimetres and inches
- Meets the single use requirement to avoid infection



**CODE**  
SMPAR

**DESCRIPTION**  
LIBERTY SKIN MARKER AND RULER  
0.5MM TIP (BOX 25)

*Capes* **CAPES**