

# Sterilisation Products for Primary Care

## Steriliser Test Bowie-Dick

Daily Air Removal - Steam Penetration

The efficiency of air removal is a critical parameter for pre-vacuum sterilisers. A steam penetration test should be performed in an empty steriliser at the beginning of every day to verify and document the efficiency of the vacuum pumps and to indicate possible leaks and/or poor steam quality.

- Get more sensitive reactions with possible air leaks with 6 test pads to indicate pass or failure
- Brass test tube enables life of 1000 cycles (with correct replacement of wear parts)
- Validated according to test pack DIN EN285



**\$393.92**  
plus GST

Order Code: F15.100.060L

Description: Bowie Dick Starter Kit with 100 Strips

Refill packs F15.100.050L Strips (pack 250) available

## Liberty Ultrasonic Bath

The affordable way to improve cleaning quality!

The unique and compact Liberty Ultrasonic Bath provides effective and reliable cleaning in a convenient table-top design. For a surprisingly affordable price! It can thoroughly clean dirt and esoposits from articles and effectively displaces contaminants. The Liberty Ultrasonic Bath is ideal for surgeries, jewellers, dentists, laboratories and domestic use.

- Featuring a stainless steel tank SUS304 with a capacity of 1L and comes complete with a basket and lid
- The unit weighs 1.6kg and measures 210x190x165mm

**\$286.00**  
plus GST



Order Code: UBL

Description: Liberty Ultrasonic Bath

**IF ITS MEDICAL, WE CAN SUPPLY**

*Distributing, representing and sourcing medical items at the best prices*

**CAPE**  
**S**

Visit **Capes NEW**  
website at

**0800 18 19 19**  
[www.capesmedical.co.nz](http://www.capesmedical.co.nz)